

BE 521: Homework 1

Exploring Neural Signals

Spring 2021

33 points

Due: Tuesday 2/2/2021 10 PM

Objective: Working with the IEEG Portal to explore different Neural signals

Shubhankar Patankar

1 Seizure Activity (16 pts)

The dataset I521_A0001.D002 contains an example of human intracranial EEG (iEEG) data displaying seizure activity. It is recorded from a single channel (2 electrode contacts) implanted in the hippocampus of a patient with temporal lobe epilepsy being evaluated for surgery. In these patients, brain tissue where seizures are seen is often resected. You will do multiple comparisons with this iEEG data and the unit activity that you worked with in Homework 0 (I521_A0001.D001). You will have to refer to that homework and/or dataset for these questions.

1. Retrieve the dataset in MATLAB using the IEEGToolbox and generate a *session* variable as before (No need to report the output this time). What is the sampling rate of this data? What is the maximum frequency of the signal content that we can resolve? (2 pts)

```
cd('/Users/sppatankar/Developer/BE-521')
addpath(genpath('ieeg-matlab-1.14.49'))
addpath(genpath('Homework_0'))
% password_file = IEEGSession.createPwdFile('spatank', '***');
session_2 = IEEGSession('I521_A0001.D002', 'spatank', 'spa-ieeglogin.bin');
sampling_rate_2 = session_2.data.sampleRate % Hz
max_frequency = sampling_rate_2/2
```

```
IEEGSETUP: Found log4j on Java classpath.
URL: https://www.ieeg.org/services
Client user: spatank
Client password: ****

sampling_rate_2 =

    200

max_frequency =

    100
```

2. How does the duration of this recording compare with the recording from HW0 (I521_A0001_D001)? (2 pts)

```
session_1 = IEEGSession('I521_A0001_D001', 'spatank', 'spa_ieeglogin.bin');
durationInUsec_1 = session_1.data(1).rawChannels(1).get_tsdetails.getDuration;
durationInSec_1 = durationInUsec_1/1e6;

durationInUsec_2 = session_2.data(1).rawChannels(1).get_tsdetails.getDuration;
durationInSec_2 = durationInUsec_2/1e6;

compare_durations = durationInSec_2/durationInSec_1
```


```
Warning: Objects of edu/upenn/cis/db/mefview/services/TimeSeriesDetails
class exist - not clearing java
Warning: Objects of edu/upenn/cis/db/mefview/services/TimeSeriesInterface
class exist - not clearing java
IEEGSETUP: Found log4j on Java classpath.
URL: https://www.ieeg.org/services
Client user: spatank
Client password: ****

compare_durations =

    64.4995
```

Recording from (I521_A0001_D002) is almost 64.5 times longer than the recording from (I521_A0001_D001).

3. Using the time-series visualization functionality of the IEEG Portal, provide a screenshot of the first 500 ms of data from this recording. (2 pts)



4. Compare the activity in this sample with the data from HW0. What differences do you notice in the amplitude and frequency characteristics? (2 pts)

The time series in (I521_A0001_D002) has a significantly higher amplitude than the one in (I521_A0001_D001). (I521_A0001_D001) has higher frequency oscillations.

5. The unit activity sample in (I521_A0001_D001) was high-pass filtered to remove low-frequency content. Assume that the seizure activity in (I521_A0001_D002) has not been high-pass filtered. Given that

the power of a frequency band scales roughly as $1/f$, how might these differences in preprocessing contribute to the differences you noted in the previous question? (There is no need to get into specific calculations here. We just want general ideas.) (3 pts)

Based on the $1/f$ scaling factor, higher frequencies have lower power. Since (I521_A0001_D001) was high-pass filtered leaving only its high frequency components, it has a generally low power. On the other hand, (I521_A0001_D002) has its lower frequency components intact implying that it has high power. Acknowledging that power and amplitude are not the same but closely related, based on their respective frequency contents, it makes sense that (I521_A0001_D001) has a lower amplitude on average than (I521_A0001_D002).

6. Two common methods of human iEEG are known as electrocorticography (ECoG) and stereoelectroencephalography (SEEG). For either of these paradigms (please indicate which you choose), find and report at least two of the following electrode characteristics: shape, material, size. Please note that exact numbers aren't required, and please cite any sources used. (3 pts)

ECoG

- Material: Platinum-Iridium
- Size: 4 mm diameter
- Source: <https://www.sciencedirect.com/science/article/pii/S105381191200403X>

7. What is a local field potential? How might the characteristics of human iEEG electrodes cause them to record local field potentials as opposed to multiunit activity, which was the signal featured in HW0 as recorded from 40 micron Pt-Ir microwire electrodes? (2 pts)

A local field potential is the transient extracellular activity recorded by an electrode from its immediate vicinity. iEEG electrodes have a non-negligible surface area that brings them in direct contact with the cortical surface resulting in recordings of local field potentials. By contrast, multi-unit activity is recorded using micro-electrodes that have a significantly smaller surface area giving them high spatial specificity.

2 Evoked Potentials (17 pts)

The data in I521_A0001_D003 contains an example of a very common type of experiment and neuronal signal, the evoked potential (EP). The data show the response of the whisker barrel cortex region of rat brain to an air puff stimulation of the whiskers. The `stim` channel shows the stimulation pattern, where the falling edge of the stimulus indicates the start of the air puff, and the rising edge indicates the end. The `ep` channel shows the corresponding evoked potential. Once again, play around with the data on the IEEG Portal, in particular paying attention to the effects of stimulation on EPs. You should observe the data with window widths of 60 secs as well as 1 sec. Again, be sure to explore the signal gain to get a more accurate picture. Finally, get a sense for how long the trials are (a constant duration) and how long the entire set of stimuli and responses are.

1. Based on your observations, should we use all of the data or omit some of it? (There's no right answer, here, just make your case either way in a few sentences.) (2 pts)

In my view, using all of the data seems reasonable. No section of the data is significantly better/worse than its other sections to justify keeping/omitting it.

2. Retrieve the `ep` and `stim` channel data in MATLAB. What is the average latency (in ms) of the peak response to the stimulus onset over all trials? (Assume stimuli occurs at exactly 1 second intervals)(3 pts)

```

session_3 = IEEGSession('I521_A0001.D003', 'spatank', 'spa_ieeglogin.bin');

sampling_rate_3 = session_3.data.sampleRate;

start_time_3 = 0;
end_time_3 = session_3.data.rawChannels(1).get_tsdetails.getEndTime/1e6; % s

ep_nr = ceil(end_time_3*sampling_rate_3);
ep = session_3.data.getvalues(1:ep_nr, 1);

stim_nr = ceil(end_time_3*sampling_rate_3);
stim = session_3.data.getvalues(1:stim_nr, 2);

window_size = 1 * sampling_rate_3; % 1 s windows
ep_windowed = reshape(ep, window_size, []);
[~, max_time_idx] = max(ep_windowed);
max_time = (max_time_idx./sampling_rate_3);
latency_ms = mean(max_time) * 1000 % ms

```

```

Warning: Objects of edu/upenn/cis/db/mefview/services/TimeSeriesDetails
class exist - not clearing java
Warning: Objects of edu/upenn/cis/db/mefview/services/TimeSeriesInterface
class exist - not clearing java
IEEGSETUP: Found log4j on Java classpath.
URL: https://www.ieeg.org/services
Client user: spatank
Client password: ****

latency_ms =

    162.1508


```

3. In neuroscience, we often need to isolate a small neural signal buried under an appreciable amount of noise. One technique to accomplish this is called the spike triggered average, sometimes called signal averaging. This technique assumes that the neural response to a repetitive stimulus is constant (or nearly so), while the noise fluctuates from trial to trial - therefore averaging the evoked response over many trials will isolate the signal and average out the noise. Construct a spike triggered average plot for the data in I521_A0001.D003. Plot the average EP in red. Using the commands `hold on` and `hold off` as well as `errorbar` and `plot`, overlay error bars at each time point on the plot to indicate the standard deviation of the responses at any given time point. Plot the standard deviation error bars in gray (RGB value: [0.7 0.7 0.7]). Make sure to give a proper legend along with your labels. (4 pts)

```

spike_trig_avg = mean(ep_windowed, 2);
time_window = linspace(0, 1, length(spike_trig_avg));
figure;
hold on
errorbar(time_window, spike_trig_avg, std(ep_windowed, [], 2), ...
    'Color', [0.7, 0.7, 0.7]);
plot(time_window, spike_trig_avg, 'r', 'LineWidth', 2)
legend('\sigma', 'STA', 'Location', 'NorthEast');
hold off
xlabel('Time from Stimulation (s)');
ylabel('\muV');
title('Spike Triggered Average');

```



4. (a) We often want to get a sense for the amplitude of the noise in a single trial. Propose a method to do this (there are a few reasonably simple methods, so no need to get too complicated). Note: do not assume that the signal averaged EP is the “true” signal and just subtract it from that of each trial, because whatever method you propose should be able to work on the signal from a single trial or from the average of the trials. (4 pts)

One approach to smoothing involves taking local averages of the signal in order to reduce the effects of noise. In order to compute the amplitude of the noise, we can remove the original signal from its smoothed version and compute the magnitude of the resulting difference.

- (b) Show with a few of the EPs (plots and/or otherwise) that your method gives reasonable results. (1 pt)

```
curr_start_time = 0;
time_window = linspace(curr_start_time, curr_start_time + 1, window_size);
noisy_signal = ep_windowed(:, curr_start_time + 1);
smoothed_signal = smoothdata(noisy_signal, 'movmean', 50);
noise = noisy_signal - smoothed_signal;
mean_noise_amplitude_t1 = mean(abs(noise))

figure;
hold on
plot(time_window, noise, 'Color', [0.7, 0.7, 0.7], 'LineWidth', 2)
plot(time_window, noisy_signal, 'b', 'LineWidth', 2)
plot(time_window, smoothed_signal, 'r', 'LineWidth', 2)
legend('Noise', 'Noisy Signal', 'Smoothed Signal', 'Location', 'NorthEast');
hold off
xlabel('Time (s)');
ylabel('\muV');
title('Smoothed Signal for Trial 1');

curr_start_time = 10;
time_window = linspace(curr_start_time, curr_start_time + 1, window_size);
noisy_signal = ep_windowed(:, curr_start_time + 1);
```

```

smoothed_signal = smoothdata(noisy_signal, 'movmean', 50);
noise = noisy_signal - smoothed_signal;
mean_noise_amplitude_t11 = mean(abs(noise))

figure;
hold on
plot(time_window, noise, 'Color', [0.7, 0.7, 0.7], 'LineWidth', 2)
plot(time_window, noisy_signal, 'b', 'LineWidth', 2)
plot(time_window, smoothed_signal, 'r', 'LineWidth', 2)
legend('Noise', 'Noisy Signal', 'Smoothed Signal', 'Location', 'NorthEast');
hold off
xlabel('Time (s)');
ylabel('\muV');
title('Smoothed Signal for Trial 11');


```


```

mean_noise_amplitude_t1 =
    414.0738

mean_noise_amplitude_t11 =
    416.4319

```





Consider trials 1 and 11. The red traces show the original signals. Blue traces are their denoised versions. For both figures, the denoised traces appear smoother. The noise traces mark the difference between the original and the smoothed signals, and in both instances seem to resemble white noise.

- (c) i. Apply your method on each individual trial and report the mean noise amplitude across all trials. (1 pt)

```
smoothed_all = smoothdata(ep_windowed, 'movmean', 50);
noise_all = ep_windowed - smoothed_all;
noise_amp_all = mean(abs(noise_all));
mean_noise_amp = mean(noise_amp_all) % uV
```

```
mean_noise_amp =
404.8957
```

- ii. Apply your method on the signal averaged EP and report its noise. (1 pt)

```
smoothed_sig_avg = smoothdata(spike_trig_avg, 'movmean', 50);
noise_sig_avg = spike_trig_avg - smoothed_sig_avg;
mean_noise_amplitude_sig_avg = mean(abs(noise_sig_avg)) % uV
```

```
mean_noise_amplitude_sig_avg =
56.0538
```

- iii. Do these two values make sense? Explain. (1 pt)

The mean noise amplitude across all trials is 404.9 μV . For the signal average EP, the noise amplitude is almost a factor of 10 lower at 56.1 μV . This makes sense since the process of taking the spike triggered average implicitly denoises the data. Applying the smoothing process to this already denoised data does not take away as much of the noise as it does when applied to individual trials.